文書番号: S516-1006

改定記号:A

外科用 X 線テレビジョン装置 WHA-200

OPESCOPE PLENO DICOM Conformance Statement



医用機器事業部

DICOM 3.0 Conformance Statement for WHA-200 Pleno

0.0 Revision History:

| Revision | Date | Description |
|----------|------------|-----------------------|
| - | 2004/02/23 | New Release |
| A | 2007/12/18 | Add Modality Worklist |

1.0 Purpose:

Define the DICOM Conformance statement associated with WHA-200.

2.0 Scope:

This document describes the DICOM Conformance statement in accordance with the document DICOM PS 3.2 Conformance.

3.0 References:

DICOM PS 3.2 Conformance

DICOM PS 3.3 Information Object Definitions

DICOM PS 3.4 Service Class Specifications

DICOM PS 3.5 Data Structures and Encoding

DICOM PS 3.6 Data Dictionary

DICOM PS 3.7 Message Exchange

DICOM PS 3.8 Network Communication Support for Message Exchange

4.0 Equipment/Materials:

N/A

5.0 Responsibilities:

N/A

6.0 Definitions:

AE – Application Entity

IOD – Information Object Definition

SCU - Service Class User

SCP - Service Class Provider

SOP – Service Object Pair

UID – Unique Identifier

7.0 Instructions:

The rest of this document is written in the format specified for DICOM Conformance statements in the DICOM PS 3.2 Conformance standard document.

8.0 Introduction

This conformance statement details the WHA-200 Pleno's compliance to DICOM 3.0. It covers all service class roles that are supported by this product:

Storage Service Class (SCU) roles

Verification Service Class (SCU) roles

Basic Grayscale Print Management Class (SCU) roles

8.1 Implementation Model

DICOM capabilities of the WHA-200 Pleno include:

The WHA-200 Pleno can send images to a remote AE by initiating the DICOM C-STORE request as a SCU.

The WHA-200 Pleno supports the DICOM Verification operation as an SCU.

The WHA-200 Pleno can send images to a DICOM Print Server AE by utilizing the services of the Basic Grayscale Print Management Meta SOP Class as an SCU.

The WHA-200 Pleno can query DICOM Modality Worklist SCP systems for patient information using the Modality Worklist Management Service Class as an SCU.

-2- S516-1006A

8.1.1 Application Data Flow Diagrams See figures 8.1.1-1, 8.1.1-2, and 8.1.1-3

Figure 8.1.1-1 Storage SCU Image Submit stored on Transfer Job Remote AE Send AE Verify Send Echo Remote AE Response Commit Request Send AE image to Storage storage Commit DICOM Standard Interface

Figure 8.1.1-2 Print SCU Queue Printer AE Print images images for printing **DICOM Standard Interface**

Query
Worklist

Modality
Worklist AE

Return
Patient Info

DICOM Standard Interface

-4- S516-1006A

8.1.2 Functional definition of Aes

Send AE:

The Send AE initiates an association with a remote AE and acts as a SCU of the Storage Service Class to store images on a remote AE that acts as a SCP of the Storage Service Class. When the image transfer is completed, the send function waits for the DIMSE-C-STORE Response from the receiving AE to indicate the status of the transfer (success or fail). When the Send AE system initiates the DICOM Echo Request, it first proposes an Association with the Verification Class Presentation Context. When the DICOM Association Accept message is received, the system sends the DIMSE-C-ECHO Request message to initiate the Verification function on the receiving AE. The status of the Verification response (success or fail) is displayed.

Print AE:

The Print AE initiates an Association with a user selected remote Print AE and acts as a SCU of the Basic Grayscale Print Management Service Class. When all of the images for a particular Film Session have been transferred, the Association is closed.

If the remote printer SCP supports the Print Job service then the Print AE can monitor the status of the Print Job on the remote printer SCP.

Modality Worklist AE:

The Modality Worklist AE initiates an Association with a user selected remote Worklist AE and acts as a SCU of the Modality Worklist Management Service Class. The Modality Worklist AE sends a C-FIND request based on parameters set by the user. The user can configure the Modality Worklist AE to query for any/all modalities supported by the local system. The user can configure the Worklist to query for exams scheduled for any AE configured in the system as a Worklist SCU. One request is sent for each modality/AE title pair configured by the user.

8.1.3 Sequencing of Real World Activities

NA

8.2 AE Specifications

8.2.1 Send AE - Specification

The Send AE provides Standard Conformance to the following DICOM V3.0 SOP Classes as a SCU:

| SOP Class Name | SOP Class UID |
|---------------------------------|------------------------------|
| Verification SOP Class | 1.2.840.10008.1.1 |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 |
| X-Ray Angiographic Image Store | 1.2.840.10008.5.1.4.1.1.12.1 |
| X-Ray RF Image Store | 1.2.840.10008.5.1.4.1.1.12.2 |

8.2.1.1 Association establishment policies

8.2.1.1.1 General

The DICOM Application Context name is 1.2.840.10008.3.1.1.1

The AE Title of the Send AE is a configurable parameter. The default title is "OEM_StoreSCU".

The Send AE establishes an association whenever a transfer job comes to the top of the transfer queue.

The Send AE establishes an association whenever the user attempts to verify the DICOM connection with a remote AE.

The maximum PDU size is 30720 bytes.

8.2.1.1.2 Number of Associations

The Send AE attempts only one Association establishment at a time.

8.2.1.1.3 Asynchronous nature

The Send AE does not perform asynchronous operations.

8.2.1.1.4 Implementation Identifying Information

The Send AE provides a single Implementation Class UID which is "1.2.840.113698.7.1".

8.2.1.2 Association initiation policy

The Send AE initiates a new association for the DIMSE-C-STORE service operation for each transfer job that comes to the top of the job queue.

The Send AE initiates a new association for the DIMSE-C-ECHO service operation.

8.2.1.2.1 Transfer Image Object to a Remote AE

8.2.1.2.1.1 Associated Real-World Activity – Queue image(s) for transfer to remote AE

The associated Real-World activity is a C-Store Request initiated by the Send AE when a transfer job comes to the top of the job queue. A transfer job is created by the user selecting an image or group of images to be sent to a remote AE.

-6- S516-1006A

8.2.1.2.1.2 Proposed presentation contexts

The Send AE proposes Presentation Contexts as shown in table 8.2.1.2.1.2-1.

The receiving AE returns which Presentation Contexts it supports in the Association Accept message.

The Secondary Capture Abstract Syntax will only be used if the receiving AE does not support any of the other proposed Abstract Syntaxes. In this case, only the modules defined for the SC IOD in <u>Annex A</u> will be supported.

Table 8.2.1.2.1.2-1 Proposed Presentation Contexts for Send AE

| Presentation Context Table | | | | | | |
|-----------------------------------|------------------------------|------------------------------------|---------------------|----------|-------------|--|
| Abstract Syntax | | Transf | Role | Extended | | |
| Name | UID | Name List | UID List | | Negotiation | |
| Secondary Capture Image Store | 1.2.840.10008.5.1.4.1.1.7 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None | |
| | | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None | |
| X-Ray Angiographic Image Store | 1.2.840.10008.5.1.4.1.1.12.1 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None | |
| | | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None | |
| X-Ray RF Image Store | 1.2.840.10008.5.1.4.1.1.12.2 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None | |
| | | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None | |
| Verification Service Class | 1.2.840.10008.1.1 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None | |
| | | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None | |

8.2.1.2.1. SOP Specific Conformance

When a successful response to a C-STORE operation is received, the status display is updated to indicate that the next image in the transfer job is being transferred.

If an Association request fails or if a Failed, Refused or Warning response to a C-STORE operation is received then the currently active transfer job is aborted from the Active transfer queue and moved to the Inactive queue.

Extended negotiation is not supported.

See Annex A for a description of the IOD modules supported.

8.2.1.2.2 Send Echo Request to Remote AE

8.2.1.2.2.1 Associated Real-World Activity - Verify DICOM connection with remote AE

The associated Real-World activity is a C-Echo Request initiated by the user to determine if a remote DICOM AE is responding.

-7-

8.2.1.2.2.2 Proposed presentation contexts

The Send AE proposes a Presentation Context as shown in table 8.2.1.2.1.2-1.

8.2.1.2.2.1 SOP Specific Conformance

The Send AE provides standard conformance to the DICOM Verification Service Class as a SCU.

8.2.1.3 Association acceptance policy

The Send AE never accepts associations.

8.2.2 Print AE - Specification

The Print AE provides Standard Conformance to the following DICOM V3.0 SOP Classes as a SCU:

| SOP Class Name | SOP Class UID |
|---|------------------------|
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 |
| Basic Film Session | 1.2.840.10008.5.1.1.1 |
| Basic Film Box | 1.2.840.10008.5.1.1.2 |
| Basic Grayscale Image Box | 1.2.840.10008.5.1.1.4 |
| Printer | 1.2.840.10008.5.1.1.16 |
| Presentation LUT | 1.2.840.10008.5.1.1.23 |
| Verification Service Class | 1.2.840.10008.1.1 |
| Print Job | 1.2.840.10008.5.1.1.14 |

8.2.2.1 Association establishment policies

8.2.2.1.1 **General**

The DICOM Application Context name is 1.2.840.10008.3.1.1.1

The AE Title of the Print AE is a configurable parameter. The default title is "OEM PrintSCU".

The Print AE establishes an association whenever a local print job comes to the top of the print queue.

The Print AE establishes an association to define a Presentation LUT for a Printer SCP.

The Print AE establishes an association whenever the user attempts to verify the DICOM connection with a remote printer AE.

The maximum PDU size is 30720 bytes.

8.2.2.1.2 Number of Associations

The Print AE can have multiple associations open at a time:

- 1. One for the Basic Grayscale Print Management service.
- 2. One for the Presentation LUT service.
- 3. Multiple associations for the Verification service.

-8- S516-1006A

8.2.2.1.3 Asynchronous nature

The Print AE does not perform asynchronous operations.

8.2.2.1.4 Implementation Identifying Information

The Print AE provides a single Implementation Class UID which is "1.2.840.113698.7.1".

8.2.2.2 Association initiation policy

The Print AE initiates a new association for the Print Service Class whenever a print job reaches the top of the print queue. The Association is closed when all of the images from the print job have been sent to the Print Server.

8.2.2.2.1 Print Image

8.2.2.2.1.1 Associated Real-World Activity – Queue images for printing

The user creates a local print job by selecting individual images or a group of images to be printed. When the local print job comes to the top of the print queue an Association Request is made. Once the Print Image Association has been established, the Print AE sends a Basic Film Session N_CREATE message to the Basic Print SCP. This is followed by a Basic Film Box N_CREATE message. The Print AE then sends a Basic Grayscale Image Box N_SET message. Finally, an N_ACTION message is sent to print images at the Basic Film Box level.

8.2.2.2.1.2 Proposed presentation contexts

The Presentation Contexts proposed by the Print AE are defined in <u>table 8.2.3.2.1.2-1</u>.

Table 8.2.2.2.1.2-1 Proposed Presentation Contexts for Print AE

| Presentation Context Table | | | | | | |
|----------------------------|------------------------|--------------------|---------------------|----------|-------------|--|
| Abstract Syntax | | Transfe | Role | Extended | | |
| Name | UID | Name List UID List | | | Negotiation | |
| Basic Grayscale | 1.2.840.10008.5.1.1.9 | DICOM Implicit | 1.2.840.10008.1.2 | SCU | None | |
| Print Management | | VR Little Endian | | | | |
| (META) | | DICOM Explicit | 1.2.840.10008.1.2.1 | SCU | None | |
| | | VR Little Endian | | | | |
| Presentation LUT | 1.2.840.10008.5.1.1.23 | DICOM Implicit | 1.2.840.10008.1.2 | SCU | None | |
| | | VR Little Endian | | | | |
| | | DICOM Explicit | 1.2.840.10008.1.2.1 | SCU | None | |
| | | VR Little Endian | | | | |
| Verification Service | 1.2.840.10008.1.1 | DICOM Implicit | 1.2.840.10008.1.2 | SCU | None | |
| Class | | VR Little Endian | | | | |
| | | DICOM Explicit | 1.2.840.10008.1.2.1 | SCU | None | |
| | | VR Little Endian | | | | |
| Print Job | 1.2.840.10008.5.1.1.14 | DICOM Implicit | 1.2.840.10008.1.2 | SCU | None | |
| | | VR Little Endian | | | | |
| | | DICOM Explicit | 1.2.840.10008.1.2.1 | SCU | None | |
| | | VR Little Endian | | | | |

8.2.2.1.2.1 SOP Specific Conformance

See Annex B for a description of the attribute values for SOP Classes proposed by the Print AE.

As individual images from the local print job are transferred to the Printer SCP the status display is updated to indicate how many images have been transferred.

If the Print Job service is supported then the Print AE can monitor the remote Print Job status and the local print job will not be removed from the print queue until a Success or Failed notification is received from the Printer SCP.

If the Print Job service is not supported then the local print job is considered completed when all of the images in the job have been transferred to the Printer SCP.

Extended negotiation is not supported.

8.2.2.2.2 Define Presentation LUT

8.2.2.2.1 Associated Real-World Activity – Select a Presentation LUT for a remote printer

The user can define a Presenation LUT that can be applied to images that are sent to the same Print AE SCP. Once the Association has been established, the Print AE sends a Presentation LUT N CREATE message to the Basic Print SCP.

8.2.2.2.2 Proposed presentation contexts

The Presentation Contexts proposed by the Print AE are defined in table 8.2.3.2.1.2-1.

8.2.2.2.2.1 SOP Specific Conformance

The Print AE provides standard conformance to the DICOM Presentation LUT Service Class as a SCU.

8.2.2.2.3 Send Echo Request to Remote Printer AE

8.2.2.2.3.1 Associated Real-World Activity – Verify DICOM connection with remote AE

The associated Real-World activity is a C-Echo Request initiated by the user to determine if a remote DICOM Printer AE is responding.

8.2.2.2.3.2 Proposed presentation contexts

The Presentation Contexts proposed by the Print AE are defined in <u>table 8.2.3.2.1.2-1</u>.

8.2.2.3.2.1 SOP Specific Conformance

The Print AE provides standard conformance to the DICOM Verification Service Class as a SCU.

The status of a C-ECHO request message is displayed (SUCCESS or FAIL).

8.2.2.3 Association acceptance policy

The Print AE never accepts associations.

-10- S516-1006A

8.2.3 Modality Worklist AE - Specification

The Modality Worklist AE provides Standard Conformance to the following DICOM V3.0 SOP Classes as a SCU:

| SOP Class Name | SOP Class UID |
|------------------------|------------------------|
| Modality Worklist Find | 1.2.840.10008.5.1.4.31 |
| Verification SOP Class | 1.2.840.10008.1.1 |

8.2.3.1 Association establishment policies

8.2.3.1.1 **General**

The DICOM Application Context name is 1.2.840.10008.3.1.1.1

The AE Title of the Modality Worklist AE is a configurable parameter. The default title is "OEM WorklistSCU".

The Modality Worklist AE establishes associations under the following conditions:

- 1. When the user initiates a manual query.
- 2. Periodically, as set up in the Auto Query configuration
- 3. When the user attempts to verify the DICOM connection with a remote Worklist AE

The maximum PDU size is 30720 bytes.

8.2.3.1.2 Number of Associations

The Modality Worklist AE can have multiple Associations open at one time:

- 1. One association establishment for each SCP configured for automatic query, and one additional association if/when the user initiates a manual query.
- 2. One association to initiate a Verification Service Echo request.

8.2.3.1.3 Asynchronous nature

The Modality Worklist AE does not perform asynchronous operations.

8.2.3.1.4 Implementation Identifying Information

The Modality Worklist AE provides a single Implementation Class UID which is "1.2.840.113698.7.1".

8.2.3.2 Association initiation policy

The Modality Worklist AE initiates a new association for the Worklist Management Class for each query session. A query session is defined as a group of queries required to completely satisfy the input from the user. The Association is closed when all of the results from the query session have been received.

The Modality Worklist AE initiates a new association to verify a DICOM connection with a remote Worklist AE when the user selects the verify option for the remote AE.

8.2.3.2.1 Worklist Query Operations

The Modality Worklist AE initiates associations to perform C-FINDs and Performed Procedure Step notifications. The association is closed after an error or when the initiator requests that it be closed.

<u>8.2.3.2.1.1 Associated Real-World Activity – Query for Scheduled Procedure information</u>

Once the Worklist Query association has been established, the Modality Worklist AE sends a series of Worklist C-FIND messages to the Worklist SCP. One C-FIND message is sent for each Modality selected by the user. One C-FIND message is also sent for each AE title selected by the user. After each C-FIND message is sent, the Modality Worklist AE waits for a C-FIND response from the SCP. If the total number of records received during the active association exceeds the maximum limit set by the user, a C-CANCEL-FIND message is sent to the SCP. Response messages are read in until a C-FIND response of Success is received. After receiving the C-FIND [Success] response, the Modality Worklist AE will send a C-FIND message for the next modality/AE Title pair. This sequence continues until all modality/AE Title pairs are queried, at which time the association is closed.

8.2.3.2.1.2 Proposed presentation contexts

The Presentation Contexts proposed by the Modality Worklist AE are defined in <u>table 8.2.4.2.1.2-1</u>=

Table 8.2.3.2.1.2-1 Proposed Presentation Contexts for Modality Worklist AE

| Presentation Context Table | | | | | | |
|-------------------------------|------------------------|------------------------------------|---------------------|-----|-------------|--|
| Abs | Abstract Syntax | | Transfer Syntax | | | |
| Name | UID | Name List | Name List UID List | | Negotiation | |
| Modality Worklist Find | 1.2.840.10008.5.1.4.31 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None | |
| | | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None | |
| Verification Service Class | 1.2.840.10008.1.1 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None | |
| | | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None | |

8.2.3.2.1.2.1 SOP Specific Conformance

The Modality Worklist AE provides standard conformance to the DICOM Modality Worklist Find Service Class as a SCU.

See Annex D for a description of the attribute values for the Modality Worklist Find operation proposed by the Modality Worklist AE.

Extended negotiation is not supported=

8.2.3.3 Association acceptance policy

The Modality Worklist AE never accepts associations.

-12- S516-1006A

8.3 Network Communication Profiles

8.3.1 Supported Communication Stacks

The WHA-200 Pleno system provides DICOM V3.0 TCP/IP Network Communication Support as defined in Part 8 of the DICOM Standard (PS 3.8).

8.3.2 OSI Stack

No OSI Stack communications are provided.

8.3.3 TCP/IP Stack

The WHA-200 Pleno system supports the TCP/IP stack.

8.3.3.2 Physical media support

The WHA-200 Pleno system is indifferent to the physical medium over which TCP/IP executes.

8.3.4 Point-to-Point Stack

No Point-to-Point Stack communications are provided.

8.4 Extensions/Specializations/Privatizations

The Storage AE (SendAE) support private attributes as defined in Annex C.

8.5 Configuration

The WHA-200 Pleno system obtains its configuration information from the following files:

merge.ini - Identifies the other three configuration files.

mergecom.pro - Defines run-time parameters.

mergecom.app - Defines services on remote AEs to which connections are possible.

mergecom.srv - Service and sequence definitions.

8.5.1 AE title/presentation address mapping

The presentation address mapping is defined in the 'mergecom.app' file. The destination AE title, host name, listen port and service list for each remote AE that the WHA-200 Pleno can connect to are defined in this file. The mapping of the hostname to an IP address is defined in the 'hosts' file.

8.5.2 Configurable Parameters

The following parameters may be configured:

In the 'mergecom.app' file:

- 1. Local AE Titles
- 2. Station name
- 3. Media storage File-Set ID
- 4. Remote AEs:
 - a. AE Title
 - b. Hostname
 - c. Port number

In the 'mergecom.pro' file:

- 1. Timeouts
 - a. Wait for Association request timeout
 - b. Wait for Association reply timeout
 - c. Wait for Association release timeout
 - d. Network write timeout
 - e. Network connect timeout
 - f. Network inactivity timeout
- 2. Maximum PDU size

In the 'hosts' file:

1. IP Addresses of remote AEs

The local network address, netmask and gateway are configured via the standard Windows Network configuration utility.

8.6 Support of Extended Character Sets

The WHA-200 Pleno supports the ISO_IR 100 Character set.

-14- S516-1006A

ANNEX A – DICOM Data Elements Supported

MODULES COMMON TO SC, XA and RF IODs

| Patient Module PS3.3 section C.7.1.1 | | | | |
|--------------------------------------|-----------|------|--|--|
| Attribute Name | Tag | Type | Description | |
| Patient's name | 0010,0010 | 2 | Patient's full legal name | |
| | | | Value is loaded from a worklist, or entered by the | |
| | | | user when creating or editing a patient list using the | |
| | | | 'New' button or the 'Edit' button. | |
| Patient ID | 0010,0020 | 2 | Primary hospital ID number or code for the patient. | |
| | | | Value is loaded from a worklist, or entered by the | |
| | | | user when creating or editing a patient list using the | |
| | | | 'New' button or the 'Edit' button. | |
| Patient's birth date | 0010,0030 | 2 | Birth date of patient | |
| | | | Value is loaded from a worklist, or entered by the | |
| | | | user when creating or editing a patient list using the | |
| | | | 'New' button or the 'Edit' button. | |
| Patient's sex | 0010,0040 | 2 | Sex of patient | |
| | | | Value is loaded from a worklist, or entered by the | |
| | | | user when creating or editing a patient list using the | |
| | | | 'New' button or the 'Edit' button. | |

| | General Study Module PS3.3 section C.7.2.1 | | | | |
|----------------------------|--|------|--|--|--|
| Attribute Name | Tag | Type | Description | | |
| Study Instance UID | 0020,000D | 1 | Unique identifier for study | | |
| | | | Value is loaded from worklist or generated by | | |
| | | | system when study is created. | | |
| Study Date | 0008,0020 | 2 | Date the Study started | | |
| | | | Value is loaded from worklist '(0040,0002) | | |
| | | | Scheduled Procedure Step Start date' or generated | | |
| | | | by system. | | |
| Study Time | 0008,0030 | 2 | Time the Study started | | |
| | | | Value is loaded from worklist '(0040,0003) | | |
| | | | Scheduled Procedure Step Start time' or generated | | |
| | | | by system. | | |
| Referring Physician's name | 0008,0090 | 2 | Patient's referring physician | | |
| | | | Value is loaded from a worklist. | | |
| Study ID | 0020,0010 | 2 | User or equipment generated Study Identifier | | |
| Accession Number | 0008,0050 | 2 | A RIS generated study number | | |
| | | | Value is loaded from a worklist, or entered by the | | |
| | | | user when creating or editing a patient list using the | | |
| | | | 'New' button or the 'Edit' button. | | |
| Study Description | 0008,1030 | 3 | User defined description of the Study | | |
| Physician of Record | 0008,1048 | 3 | Physician responsible for patient care at time of | | |
| | | | Study | | |

S516-1006A -15-

| Patient Study Module PS3.3 section C.7.2.2 | | | | | |
|--|-----------|------|--------------------------------|--|--|
| Attribute Name | Tag | Type | Description | | |
| Patient's Age | 0010,1010 | 3 | Age of the patient | | |
| Patient's Size | 0010,1020 | 3 | Height in meters | | |
| Patient's Weight | 0010,1030 | 3 | Weight in kilograms | | |
| _ | | | Value is loaded from worklist. | | |
| Occupation | 0010,2180 | 3 | Occupation of the Patient. | | |

| General Series Module PS3.3 section C.7.3.1 | | | | | |
|---|-----------|------|--|--|--|
| Attribute Name | Tag | Type | Description | | |
| Modality | 0008,0060 | 1 | Type of equipment that acquired image data (XA) | | |
| | | | Value is loaded from worklist (RF or XA) | | |
| Series instance UID | 0020,000E | 1 | Unique identifier of the Series | | |
| Series number | 0020,0011 | 2 | A number that identifies this Series | | |
| Laterality | 0020,0060 | 2C | Laterality of (paired) body part examined | | |
| Series Date | 0008,0021 | 3 | Date the Series started | | |
| Series Time | 0008,0031 | 3 | Time the Series started | | |
| Performing physician's name | 0008,1050 | 3 | Name of physician administering the Series | | |
| Protocol Name | 0018,1030 | 3 | User defined description of conditions under which | | |
| | | | Series was performed | | |
| Series Description | 0008,103E | 3 | User defined description of Series | | |
| Operator's Name | 0008,1070 | 3 | Technologist(s) supporting the Series | | |
| Body Part Examined | 0018,0015 | 3 | Text description of the part of the body | | |
| | | | examined | | |
| Patient Position | 0018,5100 | 3 | Patient position descriptor relative to the | | |
| | | | Equipment | | |
| Requested Procedure ID | 0040,1001 | 3 | ID of the Requested Procedure in the Imaging | | |
| | | | Service Request | | |
| Scheduled Procedure Step ID | 0040,0009 | 3 | ID of the Scheduled Procedure Step | | |
| Performed Procedure Step ID | 0040,0253 | 3 | ID of that part of a Procedure that has been carried | | |
| | | | out within this step | | |

| General Equipment Module PS3.3 section C.7.5.1 | | | | |
|--|-----------|------|--|--|
| Attribute Name | Tag | Type | Description | |
| Manufacturer | 0008,0070 | 2 | Manufacturer of equipment that produced images | |
| Institution name | 0008,0080 | 3 | Institution where equipment that produced images is located | |
| Institution Address | 0008,0081 | 3 | Mailing address of the institution where the equipment is located that produced the digital images | |
| Station name | 0008,1010 | 3 | User defined name identifying the machine that produced the images | |
| Manufacturer's model name | 0008,1090 | 3 | Manufacturer's model number of the equipment that produced the images | |
| Device Serial Number | 0018,1000 | 3 | Manufacturer's serial number of the equipment that produced the digital images | |
| Software version | 0018,1020 | 3 | Manufacturer's designation of software version of equipment that produced images | |

| General Image Module | | | PS3.3 section C.7.6.1 |
|-------------------------|-----------|------|--|
| Attribute Name | Tag | Type | Description |
| Image (instance) number | 0020,0013 | 2 | A number that identifies the image |
| Patient Orientation | 0020,0020 | 2C | Patient direction of the rows and columns of the |
| | | | image |
| Image (content) date | 0008,0023 | 2C | Date the image pixel data creation started |
| Image (content) time | 0008,0033 | 2C | Time the image pixel data creation started |
| Image type | 0008,0008 | 3 | See IOD specific Image Module |
| Acquisition Number | 0020,0012 | 3 | A number identifying the single continuous |
| | | | gathering of data over a period of time which |
| | | | resulted in this image |
| Acquisition Date | 0008,0022 | 3 | The date the acquisition of data that resulted in this |
| | | | image started |
| Acquisition Time | 0008,0032 | 3 | The time the acquisition of data that resulted in this |
| | | | image started |
| Images in Acquisition | 0020,1002 | 3 | Number of images that resulted from this |
| | | | acquisition of data |
| Image comments | 0020,4000 | 3 | User defined comments about image |
| Lossy Image Compression | 0028,2110 | 3 | Specifies whether an image has undergone lossy |
| | | | compression |

| | Image Pixel Module | | PS3.3 section C.7.6.3 |
|----------------------------|--------------------|------|---|
| Attribute Name | Tag | Type | Description |
| Samples per pixel | 0028,0002 | 1 | Number of samples (planes) in this image (1) |
| Photometric interpretation | 0028,0004 | 1 | Specifies the intended interpretation of the pixel data (MONOCHROME2) |
| Rows | 0028,0010 | 1 | Number of rows in image (512 or 1024) |
| Columns | 0028,0011 | 1 | Number of columns in image (512 or 1024) |
| Bits allocated | 0028,0100 | 1 | See IOD Image Module |
| Bits stored | 0028,0101 | 1 | See IOD Image Module |
| High bit | 0028,0102 | 1 | See IOD Image Module |
| Pixel representation | 0028,0103 | 1 | See IOD Image Module |
| Pixel data | 7FE0,0010 | 1 | Data stream of pixel samples which comprise the image |

| Modality | LUT Modul | e (Optio | onal) PS3.3 section C.11.1 |
|-----------------------|------------------|----------|---|
| Attribute Name | Tag | Type | Description |
| Modality LUT Sequence | 0028,3000 | 1C | Sequence of Modality LUTs (Not present if |
| | | | Rescale Intercept (0028,1052) is present) |
| > LUT Descriptor | 0028,3002 | 1C | Format of LUT Data in Sequence |
| > LUT Explanation | 0028,3003 | 3 | Free Form Text |
| > LUT Type | 0028,3004 | 1C | Specifies output values of this Modality LUT |
| > LUT Data | 0028,3006 | 1C | LUT Data (Mapping of pixel value to pixel |
| | | | intensity) |
| Rescale Intercept | 0028,1052 | 1C | Required if Modality LUT sequence is not present. |
| Rescale Slope | 0028,1053 | 1C | Required if Rescale Intercept is present. |
| Rescale Type | 0028,1054 | 1C | Required if Rescale Intercept is present. |

S516-1006A -17-

| VOI LUT Module (Optional) PS3.3 section C.11.2 | | | | |
|--|-----------|------|--|--|
| Attribute Name | Tag | Type | Description | |
| Window center | 0028,1050 | 3 | Window center for display. (512) | |
| Window width | 0028,1051 | 1C | Window width for display. Required if Window | |
| | | | center (0028,1050) is sent. (1024) | |
| VOI LUT Sequence | 0028,3010 | 3 | Sequence of VOI LUT | |
| > LUT Descriptor | 0028,3002 | 1C | Format of LUT Data in Sequence | |
| > LUT Explanation | 0028,3003 | 3 | Free Form Text | |
| > LUT Data | 0028,3006 | 1C | LUT Data (Mapping of pixel value to pixel | |
| | | | intensity) | |

-18- S516-1006A

MODULES COMMON TO XA and RF IODs

| Contrast/Bolus Module (Conditional) PS3.3 section C.7.6.4 Required if contrast media used in this image | | | | |
|--|-----------|---|-------------------------|--|
| Attribute Name Tag Type Description | | | | |
| Contras/Bolus agent | 0018,0010 | 2 | Contrast or bolus agent | |

| CINE Module (Conditional) PS3.3 section C.7.6.5 Required if pixel data is Multi-Frame Cine data | | | | |
|--|-----------|----|--|--|
| Attribute Name Tag Type Description | | | | |
| Frame time | 0018,1063 | 1C | Nominal time (msec) per individual frame. Required if Frame Increment Pointer (0028,0009) points to Frame Time. | |
| Frame time vector | 0018,1065 | 1C | An array which contains the real time increments (msec) between frames for a Multi-frame image. Required if Frame Increment Pointer (0028,0009) points to Frame Time Vector. | |
| Cine Rate | 0018,0040 | 3 | Number of frames per second | |

| Multi-Frame Module (Conditional) PS3.3 section C.7.6.6 | | | | | |
|--|-----------|------|--|--|--|
| Required if pixel data is Multi-Frame Cine data | | | | | |
| Attribute Name | Tag | Type | Description | | |
| Number of frames | 0028,0008 | 1 | Number of frames in a Multi-frame image | | |
| Frame increment pointer | 0028,0009 | 1 | Contains the Data Element Tag of the attribute | | |
| | | | which is used as the frame increment in Multi- | | |
| | | | frame pixel data. | | |

| Mask Module (Conditional) PS3.3 section C.7.6.10 Required if image may be subtracted | | | | |
|---|-----------|----|--|--|
| Attribute Name Tag Type Description | | | | |
| Mask Subtraction Sequence | 0028,6100 | 1 | Defines a sequence which describe mask subtraction operations for a multi-frame image. | |
| >Mask Operation | 0028,6101 | 1 | Identify the type of mask operation to be performed ("AVG_SUB"). | |
| >Mask Frame Numbers | 0028,6110 | 1C | Specifies the frame numbers of the pixel data used to generate the mask. | |
| Recommended Viewing Mode | 0028,1090 | 2 | Specifies recommended viewing protocols ("SUB") | |

S516-1006A -19-

| X-Ray Image Module | | | PS3.3 section C.8.7.1 |
|------------------------------|-----------|------|---|
| Attribute Name | Tag | Type | Description |
| Frame increment pointer | 0028,0009 | 1C | Required if Multi-frame image. Contains Data |
| | | | Element Tag of the attribute which is used as the |
| | | | Frame increment in Multi-frame image pixel data |
| Image type | 0008,0008 | 1 | Image identification characteristics |
| Pixel intensity relationship | 0028,1040 | 1 | The relationship between the pixel sample values |
| | | | and the X-Ray beam intensity. |
| Samples per pixel | 0028,0002 | 1 | Number of samples (planes) in the image (1) |
| Photometric interpretation | 0028,0004 | 1 | Specifies the intended interpretation of the pixel |
| | | | data (MONOCHROME2) |
| Bits allocated | 0028,0100 | 1 | Number of bits allocated for each pixel sample (8 |
| | | | or 16) |
| Bits stored | 0028,0101 | 1 | Number of bits stored for each pixel sample (8 or |
| | | | 10) |
| High bit | 0028,0102 | 1 | Most significant bit for pixel sample data (7 or 9) |
| Pixel representation | 0028,0103 | 1 | Data representation of the pixel samples (0) |

| X-Ray Acquisition Module PS3.3 section C.8.7.2 | | | |
|--|-----------|------|---|
| Attribute Name | Tag | Type | Description |
| KVP | 0018,0060 | 2 | Peak kilo voltage output of the X-Ray generator |
| | | | used |
| Tube Current | 0018,1151 | 2C | X-Ray Tube Current in mA |
| Exposure | 0018,1152 | 2C | The product of exposure time and X-Ray tube |
| | | | current expressed in mAs. Required if either |
| | | | Exposure Time (0018,1150) or X-Ray Tube |
| | | | Current (0018,1151) are not present. |
| Radiation setting | 0018,1155 | 1 | Identify the general level of X-Ray dose exposure |
| Intensifier Size | 0018,1162 | 3 | Diameter of X-ray intensifier in mm |

| Display Shutter Module (Optional) PS3.3 section C.7.6.11 | | | | |
|--|-----------|------|---|--|
| Attribute Name | Tag | Type | Description | |
| Shutter shape | 0018,1600 | 1 | Shape of the shutter defined for display (CIRCULAR) | |
| Center of circular shutter | 0018,1610 | 1C | Required if shutter shape is CIRCULAR | |
| Radius of circular shutter | 0018,1612 | 1C | Required if shutter shape is CIRCULAR | |

| X-Ray Collimator Module (Optional) PS3.3 section C.8.7.3 | | | | | |
|--|-----------|------|---|--|--|
| Attribute Name | Tag | Type | Description | | |
| Collimator shape | 0018,1700 | 1 | Shape of collimator (RECTANGULAR or | | |
| _ | | | POLYGONAL) | | |
| Collimator left vertical edge | 0018,1702 | 1C | Required if collimator shape is RECTANGULAR | | |
| Collimator right vertical edge | 0018,1704 | 1C | Required if collimator shape is RECTANGULAR | | |
| Collimator upper horizontal | 0018,1706 | 1C | Required if collimator shape is RECTANGULAR | | |
| edge | | | | | |
| Collimator lower horizontal | 0018,1708 | 1C | Required if collimator shape is RECTANGULAR | | |
| edge | | | _ | | |
| Vertices of polygonal shutter | 0018,1720 | 1C | Required if collimator shape is POLYGONAL | | |

SECONDARY CAPTURE IOD

| J | Image Pixel Module | | PS3.3 section C.7.6.3 |
|----------------------|---------------------------|---|--|
| Attribute Name | Tag Type | | Description |
| Bits allocated | 0028,0100 | 1 | Number of bits allocated for each pixel sample |
| Bits stored | 0028,0101 | 1 | Number of bits stored for each pixel sample |
| High bit | 0028,0102 | 1 | Most significant bit for pixel sample data |
| Pixel representation | 0028,0103 | 1 | Data representation of the pixel samples (0) |

| | SC Image Mo | PS3.3 section C.8.6.2 | |
|---------------------------|-------------|-----------------------|-------------------------|
| Attribute Name | Tag Type | | Description |
| Date of secondary capture | 0018,1012 | 3 | Date image was acquired |
| Time of secondary capture | 0018,1014 | 3 | Time image was acquired |

| SOP Common Module PS3.3 section C.12.1 | | | | | |
|--|-----------|---|--|--|--|
| Attribute Name | Tag | Description | | | |
| SOP class UID | 0008,0016 | Uniquely identifies the SOP class Secondary Capture | | | |
| | | Image Storage "1.2.840.10008.5.1.4.1.1.7" | | | |
| SOP instance UID | 0008,0018 | Uniquely identifies the SOP instance | | | |

X-Ray XA IOD

| X-Ray Tal | ole Module (| onal) PS3.3 section C.8.7.4 | |
|----------------|--------------|-----------------------------|------------------------|
| Re | quired if im | reated with table motion | |
| Attribute Name | Tag | Description | |
| Table motion | 0018,1134 | 2 | Is table moving or not |

| XA Positioner Module | | | PS3.3 section C.8.7.7 |
|-----------------------------|-------------|----|---|
| Attribute Name | Tag Type | | Description |
| Distance Source to Detector | 0018,1110 | 3 | Distance in mm from source to isocenter |
| Distance Source to Patient | 0018,1111 | 3 | Distance in mm from source to detector center |
| Positioner motion | 0018,1500 | 2C | Used to describe activity of imaging device |
| Positioner primary angle | 0018,1510 2 | | Position of the X-Ray image intensifier about the |
| | | | patient from the RAO to LAO direction |
| Positioner secondary angle | 0018,1511 | 2 | Position of the X-Ray image intensifier about the |
| | | | patient from the CAU to CRA direction |

| SOP Common Module PS3.3 section C.12.1 | | | | | |
|--|-----------|--|--|--|--|
| Attribute Name | Tag | Description | | | |
| SOP class UID | 0008,0016 | Uniquely identifies the SOP class X-Ray Angiographic | | | |
| | | Image Storage "1.2.840.10008.5.1.4.1.1.12.1" | | | |
| SOP instance UID | 0008,0018 | Uniquely identifies the SOP instance | | | |

X-Ray RF IOD

| SC | P Common N | Module PS3.3 section C.12.1 |
|------------------|------------|---|
| Attribute Name | Tag | Description |
| SOP class UID | 0008,0016 | Uniquely identifies the SOP class X-Ray Radiofluoroscopic |
| | | Image Storage "1.2.840.10008.5.1.4.1.1.12.2" |
| SOP instance UID | 0008,0018 | Uniquely identifies the SOP instance |

S516-1006A -21-

ANNEX B - Print AE Attributes

| SOP Class Name | Command | Attribute Name | Valid Range | Default Value |
|-----------------------|----------|-------------------------------|-------------------------------------|------------------|
| Basic Film Session | N_CREATE | Number of Copies | 1-99 | 1 |
| | | Print Priority | HIGH, MEDIUM, LOW | |
| Basic Film Session | N_ACTION | Referenced Print Job Sequence | | None |
| Basic Film Box | N_CREATE | Image Display Format | | |
| | | Film Orientation | PORTRAIT, LANDSCAPE | PORTRAIT |
| | | Magnification Type | REPLICATE, BILINEAR | None |
| | | Min Density | Depends on Printer | None |
| | | Max Density | Depends on Printer | None |
| | | Border Density | WHITE, BLACK | BLACK |
| | | Empty Image Density | WHITE, BLACK | BLACK |
| | | Trim | YES, NO | NO |
| Basic Film Box | N_ACTION | Referenced Print Job Sequence | | None |
| Basic Grayscale Image | N_SET | Image Position | 1 - 12 | Mandatory |
| Box | | Samples Per Pixel | 1 | None |
| | | Photometric Interpretation | MONOCHROME 1, MONOCHROME 2 | None |
| | | Rows | 1024 | None |
| | | Columns | 1024 | None |
| | | Pixel Aspect Ratio | 1/1 | None |
| | | Bits Allocated | 8/16 | None |
| | | Bits Stored | 8/12 | None |
| | | High Bit | 7/11 | None |
| | | Pixel Representation | 0000 | None |
| Printer | N_GET/ | Printer Status | | |
| | N_EVENT_ | Printer Status Info | | |
| | REPORT | Printer Name | | |
| | | Manufacturer | | |
| | | Manufacturer Model Name | | |
| | | Software Version | | |

-22- S516-1006A

ANNEX C - WHA-200 Pleno Private Attributes

Table C.1 Patient level private attributes

| Attribute Name | Tag | VR | VM | Description |
|--------------------|-----------|----|-----|--------------------------------------|
| Patient UID | 1011,1000 | UI | 1 | Unique identifier for Patient record |
| Miscellaneous text | 1011,1002 | LO | 1-4 | Programmable text fields |
| Equipment ID | 1011,1004 | UL | 1 | Orion Acquisition Equipment ID |
| Acquisition type | 1011,1006 | UL | 1 | Orion Acquisition type |

Table C.2 Frame attributes

| Attribute Name | Tag | VR | VM | Description |
|----------------------------|-----------|----|----|--|
| Study/Series ID | 1021,1000 | UL | 1 | System generated ID |
| Image/Frame ID | 1021,1002 | UL | 1 | System generated ID |
| Status Flag | 1021,1004 | UL | 1 | Frame status attributes |
| Frame Instance UID | 1021,1006 | UI | 1 | Unique identifier for Frame record |
| Date | 1021,1008 | DA | 1 | Date that the frame was acquired |
| Time | 1021,100A | TM | 1 | Time the frame was acquired |
| DateTime | 1021,100C | FD | 1 | Floating point representation of Frame Date/Time |
| Horizontal pixel shift | 1021,100E | FL | 1 | Sub-pixel shift of frame (column direction) |
| Vertical pixel shift | 1021,1010 | FL | 1 | Sub-pixel shift of frame (row direction) |
| Min AIO window | 1021,1012 | UL | 1 | Auto Image Optimization minimum window value |
| Max AIO window | 1021,1014 | UL | 1 | Auto Image Optimization maximum window value |
| Avg AIO window | 1021,1016 | UL | 1 | Auto Image Optimization average window value |
| Tag Fields | 1021,1018 | UL | 1 | Acquisition system generated bit settings |
| Original Study/Series ID | 1021,101A | UL | 1 | Acquisition system generated ID |
| Original Image/Frame ID | 1021,101C | UL | 1 | Acquisition system generated ID |
| Acquisition Rate | 1021,101E | FL | 1 | Rate at which frame was acquired |
| Supplemental data sequence | 1021,1020 | SQ | 1 | Sequence of supplemental data associated with |
| | | | | frame |
| > Annotation | 1021,1022 | OB | 1 | Orion frame annotation information |
| > Graphic | 1021,1024 | OB | 1 | Orion frame graphical information |
| Positioner Angle | 1021,1026 | UL | 1 | Primary positioner angle in 1/10 degrees |
| Positioner Skew | 1021,1028 | UL | 1 | Secondary positioner angle in 1/10 degrees |

S516-1006A -23-

Table C.3 Acquisition attributes

| Attribute Name | Tag | VR | VM | Description |
|------------------------------|-----------|----|-----|--|
| Next available Study ID | 1031,1000 | UL | 1 | System generated ID |
| Next available DMF ID | 1031,1002 | UL | 1 | System generated ID |
| Study ID | 1031,1004 | UL | 1 | System generated Study ID |
| Next available Series ID | 1031,1006 | UL | 1 | System generated ID |
| Acquisition type | 1031,1008 | UL | 1 | Orion acquisition type |
| Series ID | 1031,100A | UL | 1 | System generated Series ID |
| Next available Image ID | 1031,100C | UL | 1 | System generated ID |
| Original Study/Series ID | 1031,100E | UL | 1 | System generated Study/Series ID |
| Image Type | 1031,1020 | UL | 1 | System generated image type |
| Integration level | 1031,1022 | UL | 1 | Integration level used for acquisition |
| Image Study/Series ID | 1031,1024 | UL | 1 | System generated Study/Series ID stored with |
| | | - | | image record |
| Image/Frame ID | 1031,1026 | UL | 1 | System generated Image/Frame ID stored with |
| 2 | , | | | image record |
| Image status flags | 1031,1028 | UL | 1 | System generated image status bits |
| Image edge table | 1031,102A | UL | 1 | Identifies edge table used at acquisition |
| Image landmarking | 1031,102C | UL | 1 | Identifies landmarking used at acquisition |
| Image flip H/V | 1031,102E | UL | 1 | Identifies Horizontal & Vertical flipping used at |
| | , , | | | acquisition |
| Image processing default | 1031,1030 | UL | 1 | Identifies default processing applied to image at |
| settings | | | | acquisition |
| Image AIO Average goal | 1031,1032 | UL | 1 | Auto Image Optimization average goal used |
| | | | | during acquisition |
| Image AIO Maximum goal | 1031,1034 | UL | 1 | Auto Image Optimization maximum goal used |
| | | | | during acquisition |
| Image AIO Minimum goal | 1031,1036 | UL | 1 | Auto Image Optimization minimum goal used |
| | | | | during acquisition |
| LUT control points | 1031,1038 | UL | 1-8 | Control points for display LUT |
| Original image UID | 1031,103A | UI | 1 | Original UID for image |
| Digital stepping information | 1031,103C | UL | 1 | Bit settings for digital stepping information |
| Acquisition Angle | 1031,103E | UL | 1 | Acquisition angle in 1/10 degrees |
| Acquisition Skew | 1031,1040 | UL | 1 | Acquisition skew in 1/10 degrees |
| APR value | 1031,1042 | UL | 1 | Anatomical Programmed Radiology value used at |
| | | | | acquisition |
| APR table version | 1031,1044 | UL | 1 | Anatomical Programmed Radiology table version |
| Associated acquisition ID | 1031,1046 | UL | 1 | Loop/Frame ID value of an associated image (e.g. |
| | | | | bi-plane images would reference each other) |
| Rotate degrees | 1031,1048 | UL | 1 | Degrees of rotation during acquisition |
| Patient position | 1031,104A | SH | 1 | Orion code for patient position during acquisition |
| Procedure description | 1031,104C | LO | 1 | Procedure description for Orion Study record |
| Magnification Factor | 1031,104E | US | 1 | Magnification factor of acquistion system |
| Target to image distance | 1031,1050 | FL | 1 | Distance from target to imaging plane in |
| | | | | millimeters. |
| Accumulated Dose Area | 1031,1052 | UL | 1 | Study level attribute related to dosage exposure |
| Product | | | | |
| Accumulated Dose per Area | 1031,1054 | UL | 1 | Study level attribute related to dosage exposure |
| Body Surface Area | 1031,1056 | UL | 1 | Study level attribute related to dosage exposure |
| Heart Rate | 1031,1058 | UL | 1 | Heart rate in Beats per Minute (image level |
| | | | | attribute) |

-24- S516-1006A

Table C.4 Review attributes

| Attribute Name | Tag | VR | VM | Description |
|----------------------------|-----------|----|----|--|
| Shutter type | 1041,1000 | UL | 1 | Indicate type of shutter to apply to image (Auto low, Auto med, Auto high, manual) |
| Polarity | 1041,1002 | UL | 1 | Indicate if image is displayed normal or inverted |
| Edge level | 1041,1004 | UL | 1 | Edge enhancement level to apply to image |
| Zoom level | 1041,1006 | UL | 1 | Zoom factor to apply to image (2x, 3x,) |
| Zoom x/y | 1041,1008 | UL | 1 | Center point of zoom region (x in high word, y in low word) |
| Mask Image/Frame ID | 1041,100A | UL | 1 | System generated Image/Frame ID for mask image |
| Region of Interest | 1041,100C | UL | 2 | Upper/Left and Lower/Right coordinates of ROI |
| Flip Horizontal/Vertical | 1041,100E | UL | 1 | Indicate flip to apply to image (vertical in high word, horizontal in low word) |
| Loop Begin/End frames | 1041,1020 | UL | 2 | First element is start frame, second element is end frame for loop replay |
| Supplemental data sequence | 1041,1030 | SQ | 1 | Sequence of supplemental data associated with image |
| > Annotation | 1041,1032 | OB | 1 | Orion image annotation information |
| > Graphic | 1041,1034 | OB | 1 | Orion image graphical information |
| Image description | 1041,1036 | ST | 1 | User defined description of image |

S516-1006A -25-

ANNEX D - Worklist AE Attributes

Table D.1: Type 1 Fields Requested From Provider

| DICOM Tag | Description | Field Use |
|-----------|--------------------------------------|--------------|
| 0040,0100 | Scheduled Procedure Step Sequence | Sequence |
| 0040,0001 | Scheduled Station AE Title | Match |
| 0040,0002 | Scheduled Procedure Step Start Date | Match |
| 0040,0003 | Scheduled Procedure Step Start Time | Match |
| 0008,0060 | Modality | Match |
| 0040,0007 | Scheduled Procedure Step Description | Match/Return |
| 0040,0009 | Scheduled Procedure Step ID | Match/Return |
| 0040,1001 | Requested Procedure ID | Return |
| 0020,000D | Study Instance UID | Return |
| 0010,0010 | Patient Name | Match |
| 0010,0020 | Patient ID | Match |

Table D.2: Type 2 and 3 Fields Requested From Provide≠

| DICOM Tag | Description | Field Use |
|-----------|-------------------------------------|--------------|
| 0010,0030 | Patient's Birth Date | Return |
| 0010,0040 | Patient's Sex | Return |
| 0010,1030 | Patient's Weight | Return |
| 0040,0006 | Scheduled Performing Physician Name | Match/Return |
| 0008,0050 | Accession Number | Match |
| 0032,1032 | Requesting Physician | Return |
| 0008,0090 | Referring Physician | Match/Return |

-26- S516-1006A